

Unilateral Flowmeter

990622 801X-16

REPAIR PARTS

MODEL#801-OX-15

NO.	PART NO.	PART NAME	Q'TY
1	802A12	Body	1
2	FRA40	Filter	1
3	CROWA13	Holder	1
4	CROWA10H	Inner Tube	1
5	CROWA12	Safety block	1
6	CROWA09	O-Ring	1
7	CROWA11	PC Tube(OUT)	1
8	800XA20B	Aluminum Ball	1
9	800XA22A	O-Ring	1
10	800XA65	Valve Stem	1
11	800XA66A	Nut	1
12	800XA53C	Adjusting Knob	1
13	92A05	O-Ring	1

