



M40811 – Sliding Weight Chair Scale

- Robust technique with sliding weights
- The classic ADE seat, easy to clean
- Arm- and footrests swivelling
- 2 wheels
- Easy to manoeuvre by tilting the scale

The economical chair scale solution



Technical Specifications

| | |
|------------|---------------------|
| Capacity | 200 kg |
| Graduation | 100 g |
| Dimensions | 570 x 960 x 1090 mm |
| Weight | approx. 30 kg |



M40610 – Chair Scale

- Large LCD, 25 mm
- Weight unit kg/lbs selectable
- Automatic zero setting and function test
- Tare and hold function
- BMI (Body Mass Index)
- Power supply: mains and rechargeable batteries
- Arm- and footrests swivelling
- Double brake system

Colour of seat at end user's selection



| Technical Specifications | |
|--------------------------|---------------------|
| Capacity | 250 kg |
| Graduation | 100 g |
| Dimensions | 570 x 935 x 1055 mm |
| Weight | approx. 28 kg |



Optional



4 steering wheels



Seat colours: anthracite, light grey and red

M40613 – Chair Scale

- Large LCD, 25 mm
- Weight unit kg/lbs selectable
- Automatic zero setting and function test
- Tare and hold function
- BMI (Body Mass Index)
- Power supply: mains and rechargeable batteries
- Arm- and footrests swivelling, double brake system
- Upholstered seat for 100% comfort

Comfortably seated even during weighing

| Technical Specifications | |
|--------------------------|---------------------|
| Capacity | 250 kg |
| Graduation | 100 g |
| Dimensions | 570 x 935 x 1055 mm |
| Weight | ca. 28 kg |

