



Artikel-Nr. M20610

International Range

Electronic Column Scale
Large LCD, 25 mm
Automatic zero setting
Tare and hold function
BMI (Body Mass Index) function
Power supply: mains and rechargeable batteries

Accessory: MZ10019 Measuring rod

Technical Specifications

Capacity 200 kg
Graduation 100 g
Dimensions 310 x 960 x 430 mm
Weight approx. 6 kg
Platform size 310 x 300 mm



East Coast Medical



Tel: 031 2015585
072 1525142

011 462 9800
082 3762548

ADE GmbH & Co.
Hammer Steindamm 27 – 29
22089 Hamburg / GERMANY

 +49(40)432776-0
Fax +49(40)432776-10
E-mail: info@ade-hamburg.de