

M20610 – Column Scale**M30610 – Floor Scale**

- Large LCD, 25 mm
- Automatic zero setting and function test
- Tare and hold function
- BMI (Body Mass Index) function
- Power supply: mains and rechargeable batteries
- Accessory: MZ10023 Measuring rod



More flexibility - portable option:
With carry case (32610)

Technical Specifications	M20610	M30610
Capacity	200 kg	200 kg
Graduation	100 g	100 g
Dimensions (mm)	310 x 960 x 430	310 x 55 x 300
Weight approx.	5,6 kg	4,4 kg
Platform size (mm)	312 x 300	312 x 300

M20612 – Column Scale

- Easy to read LCD 24 mm
- Unique round metal platform
- Anti slip durable mat for 100% safety
- Automatic zero setting
- Tare-, Hold- and BMI function
- Telescopic height rod attachable
- Power supply: 9 V Block battery and mains adapter inclusive



All around Health Check

Technical Specifications	
Capacity	200 kg
Graduation	100 g
Dimensions	400 x 940 x 510 mm
Weight	approx. 4,5 kg
Platform size	∅ 400 mm

M35611 – Floor Scale

- Extra large and comfortable platform
- Extra sized LCD, 65 mm
- Dial framed by elegant chrome ring



The digitalized classic design

Technical Specifications	
Capacity	200 kg
Graduation	100 g
Dimensions	320 x 64 x 435 mm
Weight	approx. 2,4 kg

M20812 – Sliding Weight Column Scale

- Reading on both sides
- Robust lever system
- Platform and column: powder coated
- Weighing head: cast iron
- Platform with non-slip plastic cover
- Measuring range height rod: 750 – 2000 mm/Graduation 5 mm
- Also double graduated (kg/lbs) available



The most requested Physician Scale
All inclusive: height rod & castors

Technical Specifications

Capacity	200 kg
Graduation	100 g
Dimensions	280 x 1470 x 560 mm
Weight	approx. 15 kg
Platform size	275 x 375 mm



M20813 – Automatic Round Dial Column Scale

- Robust and massive steel case
- Easy to read dial with chrome ring
- Height rod included
- Also double graduated (kg/lbs) available
- Telescopic height rod:
measuring range: 600 - 2100 mm,
graduation 1 mm

The classic weighing and measuring station
Just stand on and read out

Technical Specifications

Capacity	150 kg
Graduation	500 g
Dimensions	280 x 1085 x 665 mm
Weight	approx. 13,5 kg
Platform size	280 x 370 mm

