

# Medical Oxygen Regulator

REPAIR PARTS

MODEL:101RM-15LM-992HT-OX

NO.	PART NO.	PART NAME	Q'TY
1	103A01HA	Body	1
2	25A163	Knob	1
3	25A17A	Inlet Nut	1
4	103A02A	Inlet Spigot	1
5	25A22	Filter	1
6	GLCWWA67	O-Ring	1
7	103A04F	Gauge	1
8	103A05F	L.P. Pressure Gauge	1
9	800A10	Safety Body	1
10	800A11	Safety Seat	1
11	800A12	Seat Retainer	1
12	800A13	Spring	1
13	800A14	Cap	1
14	106A11B	Outlet Connection	1
15	106A21	Spring	1
16	YXA17	Seat Assembly	1
17	106A05	Nozzle	1
18	C200A15	Retainer	1
19	101A76	Diaphragm	1
20	101A74B	Diaphragm Plate	1
21	101A73A	Nut	1
22	101A75	Gasket	1
23	103A20	Spring	1
24	101A68B	Spring Button	1
25	103A11	Adjusting Screw	1
26	103A10	Bonnet	1
27	103A12	Nut	1
28	101A169B	Screw	1
29	102A03	Adjusting Knob	1
30	GCWA18	Screw	1

