

# Medical Oxygen Regulator

REPAIR PARTS

MODEL:101RM-15LM-992HT-OX

NO.	PART NO.	PART NAME	QTY
1	103A01HA	Body	1
2	800A43A	Yoke Connector	1
3	103A04F	Gauge	1
4	103A05F	L.P. Pressure Gauge	1
5	800A10	Safety Body	1
6	800A11	Safety Seat	1
7	800A12	Seat Retainer	1
8	800A13	Spring	1
9	800A14	Cap	1
10	106A11B	Outlet Connection	1
11	106A21	Spring	1
12	YXA17	Seat Assembly	1
13	106A05	Nozzle	1
14	C200A15	Retainer	1
15	101A76	Diaphragm	1
16	101A74B	Diaphragm Plate	1
17	101A73A	Nut	1
18	101A75	Gasket	1
19	103A20	Spring	1
20	101A68B	Spring Button	1
21	103A11	Adjusting Screw	1
22	103A10	Bonnet	1
23	103A12	Nut	1
24	101A169B	Screw	1
25	102A03	Adjusting Knob	1
26	GCWA18	Screw	1

