

# OXIPROBE

SpO<sub>2</sub> sensor

Bio Medical Technologies(BMT) provides Pulse Oximeter sensors and cables used for non-invasive measurement of blood oxygen saturation and pulse rate measurement.


We offer both disposable and reusable sensors in compatible with several different manufacturers of pulse oximeters and patient monitors.

Choose from adhesive disposable to reusable sensor!

BMT sensors will satisfy customer need with the most cost effective solution.



CE  
0120



BIO MEDICAL TECHNOLOGIES  
[www.bm-tech.co.kr](http://www.bm-tech.co.kr)

# Quick Reference Guide

**Nellcor® Compatible Sensors** Grey Connector

Adult Disposable Finger Sensor	BM-100
Adult/Neonatal Disposable Finger/Foot Sensor	BM-200
Adult Reusable Finger Sensor	BM-300
Paediatric Disposable Finger Sensor	BM-400
Infant Disposable Big Toe Sensor	BM-500
Reusable Multi-site Sensor	BM-600

**Criticare(CSI)® Compatible Sensors** Green Connector

Adult Disposable Finger Sensor	BM-110
Adult/Neonatal Disposable Finger/Foot Sensor	BM-210
Adult Reusable Finger Sensor	BM-310
Paediatric Disposable Finger Sensor	BM-410
Infant Disposable Big Toe Sensor	BM-510
Reusable Multi-site Sensor	BM-610

**BCI® Compatible Sensors** Black Connector

Adult Disposable Finger Sensor	BM-120
Adult/Neonatal Disposable Finger/Foot Sensor	BM-220
Adult Reusable Finger Sensor	BM-320
Paediatric Disposable Finger Sensor	BM-420
Infant Disposable Big Toe Sensor	BM-520
Reusable Multi-site Sensor	BM-620

**Datex® Compatible Sensors** Yellow Connector

Adult Disposable Finger Sensor	BM-130
Adult/Neonatal Disposable Finger/Foot Sensor	BM-230
Adult Reusable Finger Sensor	BM-330
Paediatric Disposable Finger Sensor	BM-430
Infant Disposable Big Toe Sensor	BM-530
Reusable Multi-site Sensor	BM-630

**Novamatrix® Compatible Sensors** Blue Connector

Adult Disposable Finger Sensor	BM-140
Adult/Neonatal Disposable Finger/Foot Sensor	BM-240
Adult Reusable Finger Sensor	BM-340
Paediatric Disposable Finger Sensor	BM-440
Infant Disposable Big Toe Sensor	BM-540
Reusable Multi-site Sensor	BM-640

**Nonin® Compatible Sensors** Beige Connector

Adult Disposable Finger Sensor	BM-150
Adult/Neonatal Disposable Finger/Foot Sensor	BM-250
Adult Reusable Finger Sensor	BM-350

Paediatric Disposable Finger Sensor	BM-450
Infant Disposable Big Toe Sensor	BM-550
Reusable Multi-site Sensor	BM-650

**Ohmeda® Compatible Sensors** Orange Connector

Adult Disposable Finger Sensor	BM-160
Adult/Neonatal Disposable Finger/Foot Sensor	BM-260
Adult Reusable Finger Sensor	BM-360
Paediatric Disposable Finger Sensor	BM-460
Infant Disposable Big Toe Sensor	BM-560
Reusable Multi-site Sensor	BM-660

**Nihon Kohden® Compatible Sensors** Lt. Blue Connector

Adult Disposable Finger Sensor	BM-170
Adult/Neonatal Disposable Finger/Foot Sensor	BM-270
Adult Reusable Finger Sensor	BM-370
Paediatric Disposable Finger Sensor	BM-470
Infant Disposable Big Toe Sensor	BM-570
Reusable Multi-site Sensor	BM-670

**Complete Replacement Sensors with Cable**

BM-3010	Marquette	BM-3020-R	Hewlett Packard
BM-3020-D	Hewlett Packard	BM-3030	Siemens
BM-3030-R	Siemens	BM-3040	Drager
BM-3050	Spacelabs	BM-3070	Nellcor
BM-3080	Spacelabs	BM-3090	Artema
BM-3224-L	BCI	BM-3330	Datex
BM-3424	Novamatrix	BM-3630	Ohmeda
BM-3730	Nihon Kohden		

**Replacement Extension Cables**

BM-2010	Marquette	BM-2012	Nellcor
BM-2024	Nellcor	BM-2020-R	Hewlett Packard
BM-2020-D	Hewlett Packard	BM-2020-RD	Hewlett Packard
BM-2030	Siemens	BM-2030-R	Siemens
BM-2040	Drager	BM-2050	Spacelabs
BM-2070	Nellcor	BM-2080	Spacelabs
BM-2090	Artema	BM-2124-D	CSI
BM-2124-L	CSI	BM-2224	BCI
BM-2224-L	BCI	BM-2330	Datex
BM-2424	Novamatrix	BM-2512	Nonin
BM-2524	Nonin	BM-2630	Ohmeda

**For Micro-foam tape, please indicate your order number as follows:**  
M=Micro-foam e.g. BM-100M, BM-200M, BM-400M, BM-500M

# BMT Sensor Reference

		BMT Sensor	Recommended Use	Patient Size	Preferred Application	Advantage
		<b>BM-100 Series</b>	Disposable	>30kg Adult	Index finger	* Self-adhesive application
		<b>BM-200 Series</b>	Disposable	<3kg or >40kg Neonatal/Adult	Lateral aspect of foot or Index finger	* Pre-formed for easy application * Optimal performance in motion situation
		<b>BM-400 Series</b>	Disposable	10kg-50kg Pediatric	Index finger	* Single patient use supports infection control management
		<b>BM-500 Series</b>	Disposable	3kg-20kg Infant	Big toe	* Latex free
		<b>BM-300 Series</b>	Reusable	>40kg Adult	Index finger	* Fast and easy application for spot checking and short duration continuous use
		<b>BM-600 Series</b>	Reusable	>1kg	Multisite or Ear	* Latex free


- BMT Sensors provide a cost-effective
- Compatible with major leading manufacturer of patient monitoring system
- Simplified design allows for quick set-up and easy checking
- Advanced high quality of LED provides superior performance of signal tracking
- Latex free rubber based adhesive coated tape provides good adhesion to skin but removes cleanly without causing skin trauma

## BMT SpO<sub>2</sub> Sensors, Choose the right fit

BMT OxiPROBE Sensor	Replacement	Packing
BM-100 Series	BCI, Colin, Criticare, Critikon, Dixtal, Datascope, Datex, Drager, Fukuda-Denshi, Hewlett-Packard, Invivo, Kontron, Marquette, Mennen, Nellcor, Nihon-Kohden, Nonin, Novamatrix, Ohmeda, Protocol, Spacelabs, Siemens, etc.	20pcs/box
BM-200 Series		20pcs/box
BM-300 Series		1pc/box
BM-400 Series		20pcs/box
BM-500 Series		20pcs/box
BM-600 Series		1pc/box
BM-2000 Series	Extension cables and Complete sensors are available for the above models.	1pc/box
BM-3000 Series		1pc/box
BM-6000 Series		1pc/box

**Bio Medical Technologies Co., Ltd.**  
 Shin Bldg. #298-2, 5th St. Samsun-dong,  
 Sungbuk-gu, Seoul 136-045, Korea  
 Tel. 82-2-929-0540~1 Fax. 82-2-929-0542  
 E-mail: bmtshin@kornet.net


















FIME2005







































BIO MEDICAL TECHNOLOGIES  
 www.bm-tech.co.kr



















All BMT Oxiprobe sensors are CE class IIb certified. / BMT is an ISO 9001:2000/ISO 13485 certified manufacturer.










# SpO<sub>2</sub> Cross Reference List

Manufacture	Monitor model	Monitor side connector	Extension(Adapter) cable	Disposable sensor	Reusable sensor	Complete sensor	Compatible
Aqilent (Philips)	Viridia A1, A3	 DPK1	BM-2020-R DPK1 → DB9 cable 3m	BM-100	BM-300	BM-3020-R DPK1 cable 3m	Nellcor
Alaris	Vital Check 4400	 6pin DB9	BM2124-D DB9 → DB9 cable 2.4m	BM-110	BM-310	on request	Criticare (CSI)
Aqitron	Serena MX	 7pin Redel	BM2224-L Lemo → DB9 cable 2.4m	BM-120	BM-320	BM2224-L Lemo connector cable 2.4m	BCI
	9101-OX	 7pin Hypertronics	BM2630 Hypertronics → DB9 cable 3.0m	BM-160	BM-360	BM3630 Hypertronics cable 3.0m	Ohmeda
Armstrong	AD900, 1000, 2000, 3000	 9pin DB9	BM-2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI
Artema	MM200, MM204	 Artema 12pin	n/a	n/a	n/a	n/a	Artema
	Diascope 2, Gamme 8000	 7pin DB9	BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
	9050, 9060	 10pin Nicolay	BM2090 Nicolay → DB9 cable 3m	BM-100	BM-300	BM3090 Nicolay connector cable 3m	Nellcor
Baxter	S&W	 10pin Nicolay	n/a	n/a	n/a	n/a	S&W
	ASAT	 7pin Hypertronics	n/a	n/a	n/a	n/a	Simed
	Explorer, Monarch	 7pin DB9	BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
BCI	3040, 3100, 3101, 3401, 9100, 9200, Advisor, Autocorr, Fingerprint, Mini Torr, Mini Torr-Plus	 9pin DB9	BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI
	6100	 7pin Redel	BM2224-L Lemo → DB9 cable 2.4m	BM-120	BM-320	BM2224-L Lemo connector cable 2.4m	BCI
Bird	4400	 5pin Lemo	BM2124-L Lemo → DB9 cable 2.4m	BM-110	BM-310	n/a	CSI
Biosys	BPM200, 700, Guardian	 7pin DB9	BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
Bruker	40, 443, 444	 Odatam	on request	BM-120	BM-320	on request	BCI
Burdick (Spacelabs)	100, Pro2	 9pin DB9	BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI




Manufacture	Monitor model	Monitor side connector	Extension(Adapter) cable	Disposable sensor	Reusable sensor	Complete sensor	Compatible
IVY Biomedical Systems	405P, 405T	 7pin DB9	BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
	7840, 7845, Minimon7137, Micromon, Supermon, Kolormon	 12pin Nicolay	n/a	n/a	n/a	n/a	Kontron
Kontron	Same as above	 12pin Nicolay	on request	BM-100	BM-300	on request	Nellcor
	All	 9pin DB9	BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI
Lynn Medical	Mon-Mtek-30, 33, 35, 40, 80, 85, 1200, 1201, 1300, 1401, 1406, 1411, 1413, 1501, 1503, Mod02, Mod03	 9pin DB9	BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI
Magna Medical	7000 Series, Tram A.R, Tram Modular x00SL Series	 7pin Hypertronics	BM2630 Hypertronics → DB9 cable 3.0m	BM-160	BM-360	BM3630 Hypertronics cable 3.0m	Ohmeda
Marquette (GE Medical)	Eagle, Dash, Solar, Tram Modular x50SL Series	 11pin Marquette	BM2010 Marquette → DB9 cable 3m	BM-100	BM-300	BM3010 11pin Marquette cable 3.0m	Nellcor
Marquette Hellige (GE Medical)	Eagle 1000, 3000, SMU610, 611, SMK210, SMV	 Unidentified	n/a	BM-100		n/a	Nellcor
MDE	Escort I, 100, 300A, OPT11A	 9pin DB9	BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
	Escort II, OPT30, Prism, 20600, 20401-00	 12pin Lemo	on request	BM-100	BM-300	on request	Nellcor
Med Air AB	Lifesense, Pulsesense	 9pin DB9	BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI
Mediana	YM-2002 Series	 3M	BM2070 3M → DB9 cable 3.0m	BM-100	BM-300	BM3070 3M connector cable 3.0m	Nellcor
Medlab	Nanox2	 9pin DB9	BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI
	Pox10	 7pin Redel	BM2224-L Lemo → DB9 cable 2.4m	BM-120	BM-320	BM3224-L Lemo connector cable 2.4m	BCI
Medsource	BPM 200	 7pin DB9	BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor

MEK	MP-100, 110, 400, 1000		6pin DB9	Horizon XL(Ohmeda)	on request	BM-311	on request	MEK
Mennen	Horizon XL(Nellcor), Mercury1200, MPR1300		7pin Hypertronics	BM2630 Hypertronics → DB9 cable	on request	BM-160	BM3630 Hypertronics cable 3.0m	Ohmeda
	Horizon XL(Nellcor), Mercury1200, MPR1300		12pin Lemo	on request	on request	BM-100	on request	Nellcor
	Horizon XL(Nellcor), Mercury1200, MPR1300, Envoy		10pin Lemo	on request	on request	BM-100	on request	Nellcor
Nellcor	N-100, N-200		12pin Lemo	on request	on request	BM-100	on request	Nellcor
	NPB-290/5, N-390/5, N-3000, N-6000		3M	BM2070 3M → DB9 cable 3.0m	on request	BM-100	BM3070 3M connector cable 3.0m	Nellcor
	N-595, N-550		9pin DB9	n/a	n/a	n/a	n/a	Nellcor Oximax
Moriya	All(not listed above)		7pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	on request	Nellcor
	SD Series		9pin DB9	BM2224 DB9 → DB9 cable 2.4m	on request	BM-120	on request	BCI
	AL-800PA Module; BSM-1101, 2101/2, 3101, 7100, 7200, 0LV-1100A, Lifescope		DPK1	n/a	n/a	BM-170	BM-370	BM-3730 DPK1 connector cable 3m
Nihon Kohden	AL-900P Module; BSM-1102, 9510; OGS-2000		7pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	on request	Nellcor
	BSM-4104A, 4114A		3M	BM2070 3M → DB9	on request	BM-100	BM3070 3M connector	Nellcor
Nonin	BSM-2301K, 2303K, 4103, 9510J/K		Blue DB9	n/a	n/a	BM-170	on request	Nihon Kohden
	8500, 8600, 8700, 8800		7pin DB9	BM2512, 2524 DB9 → DB9	on request	BM-150	on request	Nonin
Novamatrix	8604		6pin Lemo	on request	on request	BM-150	on request	Nonin
	510, 511, 505, 515B/C, 520, 7200, 7100, Oxypleth		6pin Hypertronics	BM2424 Hypertronics → DB9 cable 2.4m	on request	BM-140	BM3424 Hypertronics cable 2.4m	Novamatrix
	500		7pin Hypertronics	on request	on request	BM-140	on request	Novamatrix
Odam	512, 513, 710, 715		6pin DB9	on request	on request	BM-140	on request	Novamatrix
	All		Odam	on request	on request	BM-120	on request	BCI
Pace Tech	Vitalmax 800 Series		5pin Switchcraft	n/a	n/a	BM-100	n/a	Nellcor
	520		5pin Switchcraft	n/a	n/a	n/a	n/a	Pace Tech

Colin	BP306, Press-mate SAT, BP88, BP508, 808, 885		7pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	on request	Nellcor
	188, 129, 129F, 128, 128F		7pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	on request	Nellcor
	118,120, 510, 511, 556		9pin Hypertronics	on request	on request	BM-300	on request	Nellcor
Corometrics (GE Medical)	511		7pin Hypertronics	BM2630 Hypertronics → DB9 cable 3.0m	BM3630 Hypertronics cable 3.0m	BM-360	BM-360	Ohmeda
	08/16		9pin DB9	BM2224 DB9 → DB9 cable 2.4m	on request	BM-120	on request	BCI
	501, 503, 504, 507, 508, 581, 1100, Poet, PONI, Scholar		5pin Lemo	BM2124-L Lemo → DB9 cable 2.4m	n/a	BM-110	BM-310	CSI
CSI	503Spot, 503DX, 504DX, 506, 506DX(NP, Mini SpO2		6pin DB9	BM2124-D DB9 → DB9 cable 2.4m	on request	BM-110	on request	CSI
	Dinamap 8700, 8710, 8720, Dinamap Plus, Oxyshuttle I, II		Compu-shield	n/a	n/a	n/a	n/a	Critikon
	Dinamap 9700, 9710, 9720, Dinamap Plus, Compact MPS, Pro 300/400, Select		7pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	BM-300	Nellcor
Datascope	Dinamap Oxytrak, 8280, 1846SX		9pin DB9	n/a	n/a	n/a	n/a	Critikon
	Accusart, Accutorr, 3/4 SAT, Accutorr Plus, Passport,		DIN 8	n/a	n/a	n/a	n/a	Datascope
	Expert, Passport(Nellcor)		7pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	on request	Nellcor
Datex-Ohmeda	Accutorr Plus(Nellcor)		3M	BM2070 3M → DB9 cable 3.0m	BM3070 3M connector cable 3.0m	BM-100	BM-300	Nellcor
	AS3, CS3, CH6, CHZ6, CHR6, GGS2, Cardiacap, II, 5, Capnomac Ultima, Light Monitor, M-Espir, Oscar II, OSP-200, Satlite Trans, Plus, ULTSVA, ULTSVIO, ULTSIO, ULTS300		DPK1	BM2330 DPK1 → DB9 cable 3.0m	BM3330 DPK1 connector cable 3.0m	BM-130	BM-330	Datex
	Ohmeda 3700, 3710, 3770, 3775, 3800, 3900, Biox 3740, 3760, 3800, 4700 Cap, 5220RGM, Rascal II, One-Tip Plus, Turfast		7pin Hypertronics	BM2630 Hypertronics → DB9 cable 3.0m	BM3630 Hypertronics cable 3.0m	BM-160	BM-360	Ohmeda
Digicare	Digiox920, P0930, Oxitrace, Digimax		7pin Reffel	BM2224-L Lemo → DB9 cable 2.4m	BM3224-L Lemo connector cable 2.4m	BM-120	BM-320	BCI
	Digimax 5500		9pin DB9	BM2224 DB9 → DB9 cable 2.4m	on request	BM-120	on request	BCI
Dima Italia	All		9pin DB9	BM2012, 2024 DB9 → DB9	on request	BM-100	on request	Nellcor

Palco	All	5pin Switchcraft Compu-shield	n/a	n/a	n/a	n/a	Palco
Physio Control	Lifepak 12		BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
Protocol (Welch Allyn)	Propaq206-EL, 204-EL, 202-EL, 106-EL Propaq242, Propaq Encore, Propaq CS, Vital Sign 300 series	 	BM2012, 2024 DB9 → DB9 n/a	BM-100	BM-300	on request	Nellcor
RGB Medical	Omicrom FT, Omicrom premium		BM2224-L Lemo → DB9 cable 2.4m	BM-120	BM-320	BM324-L Lemo connector cable 2.4m	BCI
Schiller	CM8	  	n/a	BM-100	BM-300	n/a	Nellcor
Sensor Medics	Oxyshuttle I, II	 	n/a	BM-100	BM-300	n/a	Nellcor
Siemens Medical	Sireacut700, 630 Sirem Sp02 module		BM2030 DPK1 → DB9 cable 3m	BM-100	BM-300	BM3030 DPK1 connector cable 3m	Nellcor
	SC6000, 6002XL, 7000, 8000, 9000XL		BM2030-R Round → DB9 cable 3m	BM-100	BM-300	BM3030-R Round connector cable 3m	Nellcor
	90419, 90501, 90651A-03, 90419-01, 90651A-07		BM2424 Hypertronics → DB9 cable 2.4m	BM-140	BM-340	BM3424 Hypertronics cable 2.4m	Novamatrix
	0351-D/G, 90465, 90466, 90467, 90489, 90651-07/08, 90651A		n/a	n/a	BM3424 Hypertronics cable 2.4m	n/a	Novamatrix
SpaceLabs	90351-0/6, 90465, 90466, 90467, 90489, 90651A, IM77, PC Express		BM2050 Hypertronics → DB9 cable 3.0m	BM-100	BM-300	BM3050 Hypertronics cable 3.0m	Nellcor
	90651-08, 90496 Ultraview		BM-2080 Rectangular → DB9 cable 3.0m	BM-100	BM-300	BM-3080 Rectangular cable 3.0m	Nellcor
SiMed	S-50, S-100, S-100e		n/a	n/a	BM3424 Hypertronics cable 2.4m	n/a	Novamatrix
	90651-08, 90496 Ultraview		n/a	n/a	n/a	n/a	Novamatrix
Vitalog	All	 	BM2124-L Lemo → DB9 cable 2.4m	BM-110	BM-310	n/a	CSI
Weimann	Oxycount mini		BM2224 DB9 → DB9 cable 2.4m	BM-120	BM-320	on request	BCI

Subject to change without notice

Dixtal	DX2010 series, DX8100		BM2424 Hypertronics → DB9 cable 2.4m	BM-140	BIV-340	BM3424 Hypertronics cable 2.4m	Novamatrix
Dräger	Oxipac 2000, 2500 Vitalert 1000, Cicero, Narkomed, Vitara, PM8014, 8060	  	BM2224 DB9 → DB9 cable 2.4m BM2040 Lemo → DB9 cable 3.0m	BM-120 BM-110	BM-320 BIV-310	on request BM3040 Lemo connector cable 3.0m	BCI Nellcor
Elmed	Narkomed PSO1000	 	on request BM2224 DB9 → DB9 cable 2.4m	BM-100 BM-120	BM-300 BIV-320	on request on request	Nellcor BCI
Entri Medical	Propaq102,104,106	 	on request	BM-100	BM-300	on request	Nellcor
Fukuda Denshi	HG302 Module, DS5100E, DS-5300W		BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
Gambro	BP806		BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
Getemed	VitalGuard		BM2012, 2024 DB9 → DB9	BM-100	BM-300	on request	Nellcor
Goldway	UT4000A, PS01000 UT6000A/F	 	BM2224-L Lemo → DB9 cable 2.4m BM2012, 2024 DB9 → DB9	BM-120 BM-100	BM-320 BIV-300	BM3224-L Lemo connector cable 2.4m on request	BCI Nellcor
Healthdyne	All		BM2512, 2524 DB9 → DB9	BM-150	BIV-350	on request	Nonin
Hewlett Packard (Philips)	7835C, 7835C, 7883A/C, M1020A, M1025B, M1205A, M1350B, M1722A/B, M1732A/B, M2475B	 	BM2020-R DPK1 → DB9 cable 3m BM2020-D D shape → DB9 cable 3m	BM-100 BM-100	BM-300 BIV-300	BM3020-R DPK1 cable 3m BM3020-D D shape cable 3m	Nellcor Nellcor
Innomed	Oxycard/Innocare-T		BM2512, 2524 DB9 → DB9	BM-150	BM-350	on request	Nonin
Invacare	All	 	BM2224 DB9 → DB9 cable 2.4m	BM-120	BIV-320	on request	BCI
Invivo Research	4500, 4500Plus, Millermia, Omega1445, Omni Trak TVS, Sat Scout	 	n/a	n/a	n/a	n/a	Invivo
IVV Biomedical Systems	Same as above and 4500 Plus2 405A, 405D	 	on request on request	BM-100 BM-100	BIV-300 BIV-300	on request on request	Nellcor Nellcor