

Connecting



from
heart
to
heart



SonoTrax Series Ultrasonic Pocket Doppler can meet the routine examination requirements of obstetricians, midwives and pregnant women.

SonoTrax

Ultrasonic Pocket Doppler

Features

- Easy-to-use
- Ergonomic design
- Interchangeable probe
- Built-in audio recorder
- Large backlight LCD display
- Accurate FHR detection with clear sound
- Support 2MHz, 3MHz, 4MHz, 8MHz waterproof probes
- Automatic and manual FHR counting
- Long time continuously use (4 hours)
- Automatic power-off



SonoTrax

Ultrasonic Pocket Doppler

Support 2MHz, 3MHz,
4MHz, 8MHz Probes



Backlight LCD



Built-in Recorder



Interchangeable Waterproof Probe



Carry Bag



Specifications

Display: 45x25mm LCD display

Heart rate range: 50~210bpm

Resolution: 1bpm

Audio sampling frequency: 4KHz /8KHz

Recording length: 480s/240s

Audio output power: 0.5W

Auto shut-off: 1 minute

Size: 32 (Depth) x 85 (Width) x 138 (Height) mm

Weight: 290±5g (including one battery)

Model	Interchangeable probes	FHR display	Built-in audio recorder	Backlight	Waterproof probe	Probe auto-detection display	Battery status display	Low battery indication
SonoTrax Lite	Yes	No	No	No	Yes	No	No	Yes
SonoTrax Basic	Yes	Yes	No	No	Yes	No	No	Yes
SonoTrax Basic A	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes
SonoTrax Pro	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Optional accessories	Rechargeable Li-ion battery, charger, earphone, carry bag, 9V alkaline battery							