

Suction Liners

Product Features

1. The disposable suction liner is manufactured from polyethylene which is safe and non-toxic, meeting environmental protection requirements
2. The suction liner has an optional solidifier capsule that will solidify liquids, making the handling of the bags easier and avoids accidental spillage.
3. The suction liner has a reliable seal that can stand a 60Kg pressure.
4. Double safety design - the suction liner has a built in filter and a non-return valve to ensure no fluid can be accidentally sucked into the suction controller, or central vacuum line.
4. 1000ml, 2000ml, 3000ml liners available, with various lid options - depending on the application.



3000ml
3.0L Liner shut-off
valve Lid



3000ml
3.0L Liner without
shut-off valve Lid



2000ml
2.0L Liner shut-off
valve Lid



2000ml
2.0L Liner without
shut-off valve Lid



Product Code No.	Description
G200934	2.0L Suction Liner with Filter(Aqua lid)
G200935	2.0L Suction Linerwith Filter(Aqua lid) with Pre-loaded Solidifier
G200936	2.0L Suction Liner without Filter(White lid)
G100934	1.0 L Suction Liner with Filter(Aqua lid)
G301132	3.0 L Suction Liner with filter(Aqua Lid)
G301133	3.0L Suction Liner with filter (Aqua Lid) with Pre-loaded Solidifier
G301134	3.0 L Suction Liner without filter(White Lid)
G301135	3.0 L Suction Liner(White Lid) with Pre-loaded Solidifier
S12-03	3.0 L Canister
S12-02	2.0 L Canister
S12-01B	1.5 L Canister
S12-04	1.0 L Canister
G10235	External Shutoff Filter/ VAC-GARD