



COAGULATOR^{Plus}

Procedural Versatility

With a multitude of applications in such diverse specialties as dermatology, plastic surgery and general surgery, the Coagulator Prime lets you perform a wide variety of desiccation, fulgration and coagulation office procedures.

Easily Controllable Power Setting

Digitalized display and control enhance your skill, accuracy and control when performing delicate procedures.

Reusable Hand Piece

By using reusable hand control switch, you can set power easily, and you can change variety of electrodes.

Cost-Effective

STANDARD ACCESSORIES



■ Autoclavable handpiece



■ Attractive compact electrodes



■ Power Cord



■ Foot Switch

SPECIFICATIONS

- ◆ Input Power 110V 1A 50~60Hz
220V 0.5A 50~60Hz
- ◆ Weight 2.2kg
- ◆ Dimensions 230(L) × 187(W) × 85(H)mm

- ◆ Output Characteristics
 - FREQUENCY : 500KHz
 - MONOPOLAR : 30W (500 Ω)
 - BIPOLAR : 20W (100 Ω)

*Specifications are subject to change without notice for improvement.

UM-D30 Plus (3) 2001 6.01



**UNION MEDICAL
ENGINEERING CO., LTD**

#71, Kwangsa-Ri, Yangju-Eub, Yangju-Kun,
Kyengki-Do, Korea. 482-820
Tel:82-31-840-5906-8 Fax:82-31-845-2283
www.union-medical.com