

M50610 – Wheelchair Scale

- Large LCD, 25 mm
- Unit function kg/lbs selectable
- Automatic zero setting and function test
- Tare and hold function
- BMI (Body Mass Index) function
- Operated by mains and rechargeable batteries
- Can be stored vertically
- One ramp included



Technical Specifications

Capacity	300 kg
Graduation	100 g
Dimensions	915 x 44 x 800 mm
Weight	approx. 22 kg
Platform size:	800 x 800 mm

Easy to manoeuvre on 2 castors
Suitable for electronic wheelchairs



Options:



2nd ramp MZ30010



Secure side boundaries
MZ30012



Scale with foldable column
M50613



M50612 – Wheelchair Scale

- Large LCD, 25 mm
- Unit function kg/lbs selectable
- Automatic zero setting and function test
- Tare and hold function
- BMI (Body Mass Index) function
- Power supply:
mains and rechargeable batteries



ADE Lightweight

Technical Specifications

Capacity	200 kg
Graduation	100 g
Dimensions	740 x 40 x 900 mm
Weight	approx. 15 kg

Easy transport on 2 castors
Needs no space as to be stored vertically