

Model

YP-90

YP-90B

Standard configuration

YP-90: Main body (including the hood, temperature controller, infant bed), I.V. Pole, mattress, stand.

YP-90B: Base on YP-90, added with the phototherapy (blue light) radiant box .

Optional parts: Baby scale.



Specification

- Power supply : AC220V/ 50HZ
- Input power : 600VA
- Mode of control : Air-temperature control
- Control range : 25 ℃ ~37 ℃
- Variability of temperature : < 0.5 ℃
- Uniformity of temperature : < 0.8 ℃
- Inclination of the infant bed : ± 5 °
- Noise inside hood : < 55dB (A)
- Failure alarm : over-temp alarm, deviation alarm, sensor alarm, fan motor alarm, power failure alarm, and so on.
- Average radiation of blue light : > 500 μ W/cm² (YP-90B)
- Timer range : 0-9999hrs 59min (YP-90B)

Package

Package : each unit is packed to one carton

Size : 1085 mm × 640mm × 1110mm

Gross weight: 80 Kg (YP-90) / 90 Kg (YP-90B)