

M10618 – Baby Scale

- Large LCD, 30 mm
- Unit function kg/lbs selectable
- Tare and multiple tare function
- Memory function
- Automatic shut off
- Battery use or mains adapter (included)



Safe and comfortably extra large

Carry bag MZ10025

Technical Specifications	
Capacity	20 kg
Graduation	5 g
Dimensions	600 x 385 x 85 mm
Weight	2,6 kg
Baby tray	600 x 280 x 40 mm



M10614 – Baby/Toddler Scale

- Removable baby tray
- Large LCD, 25 mm
- Unit function kg/lbs selectable
- Tare function
- Memory and recall function
- Difference value between 2 weighings
- Power supply: Battery
- Optional with mains adapter
- Available outside EC only



Helps to control baby's feeding

Technical Specifications	
Capacity	20 kg
Graduation	10 g
Dimensions	560 x 105 x 370 mm
Weight	approx. 3,5 kg
Baby tray	560 x 105 x 275 mm

M10616 – Baby/Toddler Scale

- Removable baby tray
- Tare and automatic hold function
- Unit function kg/lbs selectable
- Difference value between 2 weighings
- Automatic shut off, Power supply: Battery



Keeps mother & baby happy – every single day

Technical Specifications		M10616
Capacity		20 kg
Graduation		10 g
Dimensions (mm) approx.	525 x 143 x 290	
Weight approx.		2,4 kg
Baby tray (mm) approx.	525 x 90 x 275	



M10615 – Baby Scale

- Large LCD, 20 mm
- Weight unit kg/lbs selectable
- Automatic zero setting
- Tare and hold function
- Power supply: mains and rechargeable batteries



Technical Specifications

Capacity	20 kg
Graduation	1 g
Dimensions	535 x 200 x 330 mm
Weight	approx. 4,1 kg
Baby tray	535 x 60 x 275 mm

The typical hospital scale

M10813 – Sliding Weight Baby Scale

- Very solid steel construction
- Robust lever system
- High precision sliding weights
- Safe and easy to clean baby tray



The new edition of a classic design
For each hospital imperative

Technical Specifications

Capacity	16 kg
Graduation	10 g
Dimensions	575 x 200 x 300 mm
Weight	approx. 5,8 kg
Baby tray	575 x 70 x 275 mm

M10812 – Sliding Weight Baby Scale

- Removable baby tray
- Extra terry cloth cover for tray
- Material baby tray stainless steel, lacquered
- Easy to read sliding weights



The most economic sliding weight baby scale

Technical Specifications

Capacity	13 kg
Graduation	10 g
Dimensions	522 x 120 x 282 mm
Weight	approx. 2,3 kg
Baby tray	522 x 65 x 275

M10811 – Automatic Round Dial Baby/Toddler Scale

- Removable baby tray
- Manual zero adjustment
- Tare function
- Also double graduated (kg/lbs) available

**Technical Specifications**

Capacity	20 kg
Graduation	up to 10 kg: 50 g up to 20 kg: 100 g
Dimensions	588 x 343 x 559 mm
Weight	approx. 4,5 kg
Baby tray	550 x 265 mm

Two indicators for double reading range

M15610 – Electronic Baby Hanging Scale

- Large LCD, 15 mm
- Weight unit kg/lbs selectable
- Automatic zero setting and switch off
- Tare function
- Low battery indication
- Power supply: 2 lithium batteries



Preferred by mid-wives world-wide

Technical Specifications	
Capacity	20 kg
Graduation	20 g
Dimensions	55 x 100 x 22 mm
Weight	approx. 200 g

M15611 – Electronic Baby Hanging Scale

- LCD display
- Tare and hold function
- 2 solid hooks, 1 carabiner
- Power supply: Battery



Classic design guarantees safety and comfort

Technical Specifications	
Capacity	20 kg
Graduation	10 g
Dimensions	80 x 150 x 45 mm
Weight	approx. 500 g

Mechanical Baby Hanging Scale

M15810 – 10 kg/50 g, **M15811** – 25 kg/100 g

- Solid metal construction
- 2 solid safety hooks
- Readability: 25 g/50 g
- With 1 pair of weighing trousers/sling
- Accessories: extra trousers/carrying bag
- Also double graduated (kg/lbs) available



The most relied on baby scale throughout the world

Technical Specifications	
Capacity	10/25 kg
Graduation	50/100 g
Dimensions	Ø 170 mm
Weight	approx. 1,5 kg

Carry bag MZ10022



Weighing trousers MZ10018

