

RT2-2016 Product Catalogue



INDEX

Category	Page
Flowmeters	2
Suction Units	3
Cardiotocograph	4
Fetal Doppler	5
Electrocardiograph	6
Vital Signs Monitors & Accessories	7-9
Patient Monitors & Accessories	10-12
ICU Modular Monitors	13
Electrosurgical & Accessories	14-15
Incubators	16



RT2 Item Code

41112501-00001

OXYGEN FLOWMETER - SINGLE

FEATURES:

- Precision engineered pressure compensated design to accurately measure flow from 1 to 15 LPM
- Easy to read flow tube manufactured from crack resistant polycarbonate.
- Chromed brass body for durability
- Stainless steel needle valve for precise control
- Metal control Knob for durability (as per RT2 specification)



Each Flow Meter supplied with a metal tail pipe, as per RT2 specification

RT2 Item Code

41112501-00002

OXYGEN FLOWMETER - DOUBLE

FEATURES:

- Precision engineered pressure compensated design to accurately measure flow from 1 to 15 LPM
- Easy to read flow tube manufactured from crack resistant polycarbonate.
- Chromed brass body for durability
- Stainless steel needle valve for precise control
- Metal control Knob for durability (as per RT2 specification)



Each Flow Meter supplied with a metal tail pipe, as per RT2 specification

PLEASE NOTE

Above Flow Meters plug directly into the wall Oxygen Outlet.
Should a Flow Meter be required to fit onto Gabler Rail, then the following item will be required: (additional cost)

RT2 Item Code

42221707-00001

GABLER SNAP LOCK INCLUDING OXYGEN HOSE AND HAYER OR DRAGER CONNECTION, REQUIRED FOR OXYGEN FLOWMETER, DOUBLE



200ml Humidifier Bottle, with crack resistant polycarbonate bottle and lid
2 psi safety valve
Re-usable

RT2 Item Code

42221707-00002

The most popular high performance model of Victoria especially suitable for wards, operating rooms, recovery rooms, ICUs, obstetrics and endoscopic applications, etc.

RT2 Item Code

42171801-00002

**SUCTION UNIT, MOBILE,
MOUNTED ON MOBILE STAND**

Victoria Suction Pump

- Makes less noise than your PC!
A special system of noise and vibration dampening for an extremely quiet and smooth operation with an incredibly low noise level of 39,4 dB
- Maintenance-free Motor
- Each pump is thoroughly tested prior to delivery
- Designed for Continuous operation (100% Duty Cycle Rated)
- Casing is made of very sturdy, unbreakable, fire retardant material with UV protection

Unit can be used for drainage, where central suction is not available. Due to the quiet operation, the unit will not disturb a patients sleep. Rated to run for weeks, without switching the pump off.



The following accessories supplied with the unit:

Re-Usable Bottles with lids - Autoclaveable *
Silicone tubes set 8 x 14 mm
Tube connectors
Probe connector
Aspiration probe.

RT2 Item Code

42171803-00003

**SUCTION UNIT, WALL-MOUNT, A COMPLETE SINGLE
SUCTION UNIT WITH AN ACCESSORY RAIL
MOUNTING CLIP.**



Part Number: D2-002A-603-300

Low Vacuum Regulator (Max 300 mmHg)

Part Number: D2-002A-603-760

High Vacuum Regulator (Max 760 mmHg)

FEATURES:

- Dual Polycarbonate Canister System
- Rail Mount Fittings
- 1.5 Meter Suction hose to Wall Point
- Supplied with disposable liners
- 2 x 2 Ltr canister supplied as standard

RT2 Item Code

42181701-00005

CARDIOTOCOGRAPH, TWINS MONITORING

With its elegant design, bright colour screen, and enhanced functionality, the F3 Fetal monitor will address the needs of obstetric departments in both clinics and hospitals. F3 is amazingly lightweight and portable, which makes it suitable for outpatient service. With the latest technologies from EDAN, F3 offers an extensive set of monitoring parameters such as FHR, TOCO, DECG, IUP and fetal movement

F3

Fetal Monitor



Multiple Display Modes

Switching among three display modes, user may choose the most suitable one for clinical use.



Supplied with mobile cart

Features

- Compact, portable and lightweight
- Long life lithium-ion battery for portability
- 5.6" foldable colour TFT screen
- Large numeric and waveform display
- Simple workflow and Advanced FHR Detection Technology
- Signals Overlap Verification to differentiate twins/maternal FHR
- FHR signal quality indicator helps optimize probe position
- Built-in thermal printer
- Print on 150/152 mm width recording paper
- 1/2/3 cm/min real-time printing speed
- 15 mm/sec fast printing speed for history traces

Powerful Data Management

- 12 hours built-up memory for seamless monitoring

Accessories Supplied with unit

Ultrasound transducers x 2
Toco transducer x 1
Transducer belts x 3
Marker Switch
Rechargeable batteries and mains cable
Consumables: CTG Paper x 3 Packs, Ultrasound Gel (250ml) x 1



EDAN
A world of potential



FETAL DOPPLER

RT2 Item Code

42201709-00001

MONITOR, DOPPLER, FETAL, HEART RATE

SonoTrax Series Ultrasonic Pocket Dopplers' superior sensitivity, high-fidelity sound, recording option and optional rechargeable battery ensure the doppler family to meet the routine examination requirements of obstetricians and midwives.



Features

- Ultrasonic Pocket Doppler with Audio Recording and Playback.
- Easy to use.
- Compact design.
- Back-lit LCD Display.
- 2MHz Interchangeable waterproof probe.
- Accurate FHR detection with high-fidelity sound effect.
- Audio recorder for easy playback at any time. (240s)
- Supports 2MHz, 3MHz, 4MHz, 5MHz, 8MHz waterproof probes.
- Auto, Average and Manual modes.
- Re-Chargeable Battery lasts 4 hours of continuous work.
- Battery charge indicator.
- Automatic power-off function.
- Supplied with case and probe

8 Mhz probe can be used for vascular work (optional)



Compact Battery Charger

Backlit LCD

Built-in Recorder

Interchangeable Waterproof Probe



Supplied with 2MHz Probe, rechargeable batteries, charger and carry bag as per RT2 specifications.





RT2 Item Code

42181701-00006

SE-1200 Express

12 Channel ECG Machine

ELECTROCARDIOGRAPH (12 CHANNEL) WITH TROLLEY

The SE-1200 Express Electrocardiograph is a high-end and easy-to-operate system that provides basic 12-lead resting ECG analysis, reliably and accurately. A complete suite of advanced ECG analysis programs make the SE-1200 versatile enough to handle a wide range of cardiac diagnostic tests.



The following accessories supplied with the unit

Rolling stand with drawers
ECG cable (with set of adaptor clips for disposable electrodes)
Suction Cups - set of 6
Limb Clamps - set of 4
Pack of ECG Paper
Pack of disposable electrodes (50 electrodes)
Mains Cable
Operator's manual

FEATURES

- 8" High Resolution Colour Touch Screen
- Alphanumeric keyboard and one touch operation
- Auto Analyzing and Diagnostic Software (Interpretation)
- Real-time waveform freezing
- 120 seconds ECG waveform reviewable
- Heart rate variability (HRV) analysis
- Built-in rechargeable Li-ion battery
- 200 ECGs in internal memory
- Record 12 ECG traces on A4 size or standard US Letter-size paper
- Pre-sample printing, Trigger-sample printing for arrhythmia
- Waveform freeze, review and edit
- Diagnose results and preview before printing
- Battery capacity: 4Hr monitoring or 200 ECG's
- AC Supply: 100V-240V, 50/60HZ
- Data transfer via Ethernet and RS232
- Accurate ECG Interpretation Program
- ST Segment analysis



Unit is supplied with an Elite Rolling Stand:

- Stable Cast Iron Base.
- Locking Castors
- 2 x Drawers for paper and accessories
- ECG Machine secured to trolley, to avoid unauthorised removal of ECG Machine.



Built in LAN Card
Save files in PDF Format for Transmission and Storage





EDAN
A world of potential



M3A



RT2 Item Code

42181701-00004

Vital Signs Monitor
NIBP / SPO2 / Pulse Rate

EDAN M3A Series of Vital Signs Monitor from EDAN Instruments, provides essential measurements for patient's SpO2, NIBP and Pulse Rate. It offers a high cost-performance ratio to give practitioners throughout the hospital, most of the information they need at the patient's side.



- 3.5" colour TFT LCD and clear LED numeric display
- NIBP, SpO2, Pulse Rate Measurements
- Battery working time up to 21 hours
- Spot-check and continual working modes
- Powerful storage capacity
- Suitable for adult, pediatric and neonatal patients

The following accessories are supplied with the unit:

Adult NIBP cuff
Air hose
SpO ₂ Finger sensor for adults with cable.
Battery



RT2 Item Code

42272501-00001

Vital Signs Monitor
NIBP / SPO2 / Pulse Rate on Trolley



- 3.5" colour TFT LCD and clear LED numeric display
- NIBP, SpO2, Pulse Rate Measurements
- Battery working time up to 21 hours
- Spot-check and continual working modes
- Powerful storage capacity
- Suitable for adult, pediatric and neonatal patients

The following accessories are supplied with the unit:

Adult NIBP cuff
Air hose
SpO ₂ Finger sensor for adults with cable.
Battery



RT2 Item Code

42272501-00002

Vital Signs Monitor
NIBP / Pulse Rate



- 3.5" colour TFT LCD and clear LED numeric display
- NIBP, Pulse Rate Measurements
- Battery working time up to 21 hours
- Spot-check and continual working modes
- Powerful storage capacity
- Suitable for adult, pediatric and neonatal patients

The following accessories are supplied with the unit:

Adult NIBP cuff
Air hose
Battery



For comprehensive brochures, including full specifications, please see our website www.eastcoastmedical.co.za

Western Cape: 021 905 7721, Eastern Cape: 043 736 3111, KZN: 031 201 5585, Gauteng: 012 661 3513



RT2 Item Code

42272501-00003

Vital Signs Monitor

NIBP / Pulse Rate with trolley

EDAN M3A Series of Vital Signs Monitor from EDAN Instruments, provides essential measurements for patient's SpO₂, NIBP and Pulse Rate. It offers a high cost-performance ratio to give practitioners throughout the hospital, most of the information they need at the patient's side.



- 3.5" colour TFT LCD and clear LED numeric display
- NIBP, Pulse Rate Measurements
- Battery working time up to 21 hours
- Spot-check and continual working modes
- Powerful storage capacity
- Suitable for adult, pediatric and neonatal patients

The following accessories are supplied with the unit:

Adult NIBP cuff
Air hose
Battery



RT2 Item Code

42181801-00001

Vital Signs Monitor

SPO₂ / Pulse Rate

- 3.5" colour TFT LCD and clear LED numeric display
- SpO₂, Pulse Rate Measurements
- Battery working time up to 21 hours
- Spot-check and continual working modes
- Powerful storage capacity
- Suitable for adult, pediatric and neonatal patients



The following accessories are supplied with the unit:

SpO ₂ Finger sensor for adults with cable.
Battery



RT2 Item Code

42181801-00002

Vital Signs Monitor

SPO₂ / Pulse Rate on Trolley

- 3.5" colour TFT LCD and clear LED numeric display
- SpO₂, Pulse Rate Measurements
- Battery working time up to 21 hours
- Spot-check and continual working modes
- Powerful storage capacity
- Suitable for adult, pediatric and neonatal patients



The following accessories are supplied with the unit:


SpO ₂ Finger sensor for adults with cable.
Battery





Vital Signs Monitor Accessories Available on RT-2016 Contract

PIC	RT2 Item Code	Accessories
1	42181716-00048	CUFF, NIBP, ADULT, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
2	42181716-00049	HOSE, AIR, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
3	42181716-00050	SENSOR, FINGER, SPO2, ADULT, WITH CABLE, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
4	42181716-00051	BATTERY, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
5	42181716-00052	CUFF, ADULT, OBESE, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
6	42181716-00053	CUFF, PAEDIATRIC, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
7	42181716-00054	CUFF, NEONATAL, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
8	42181716-00056	SENSOR, SPO2, FINGER, PAEDIATRIC, WITH CABLE, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
9	42181716-00057	SENSOR, EARLOBE, ADULT, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
10	42181716-00058	SENSOR, EARLOBE, PAEDIATRIC, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2
11	42181716-00059	SENSOR, DISPOSABLE, ADULT, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2, PACK OF TEN,
12	42181716-00060	SENSOR, DISPOSABLE, PAEDIATRIC, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2, PACK OF TEN,
13	42181716-00061	SENSOR, DISPOSABLE, NEONATAL, REQUIRED FOR MONITOR, DESKTOP, NIBP AND SPO2, PACK OF TEN,





The new EDAN iM50 patient monitor provides a perfect transportation solution. The multi-functional all-in-one parameter module and new application platform greatly enhance patient monitoring.

ECG, SPO₂, NIBP, RESPIRATION, TEMPERATURE (IBP and CO₂ optional)

RT2 Item Code

42181701-00001

Offer 1



Li-Ion Battery



Portability



Bed Rail Hook

iM50
Patient Monitor

- 8.4" colour TFT display with maximum 11 waveforms
- Pacemaker detection, electrosurgical interference proof
- Defibrillation protection and defibrillation synchronization
- Full Touchscreen
- OxyCRG function for neonates
- Large font display
- Arrhythmia and S-T segment analysis
- Built-in rechargeable Li-ion battery
- USB and serial ports enables abundant future upgrades
- VGA output and analog output
- Nurse Call function and bi-directional communication with MFM-CMS central station

Standard Accessories supplied with Basic Monitor.	Quantity
ECG Cable, 3 lead, snap, Defib, IEC, 3.5m, reusable	1
ECG Electrodes, adult, disposable, 30 pieces	30pcs
SPO ₂ Finger Sensor, adult, 2.5m, reusable (can be exchanged for child or Paed Probe)	1
NIBP Cuff, Adult	1
NIBP Cuff, Large Adult	1
NIBP Cuff, Child	1
NIBP Tube, 3m	1
Temperature Probe Skin, adult, 3m, reusable	1
Power Cord (3 Mtr)	1
Rechargeable Lithium-Ion Battery (14.8V, 4200mAh)	1
User Manual	1
Ground Cable	1

Full wall mounting and rolling stand options available

RT2 Item Code

42181701-00001

Offer 2



iM60
Patient Monitor

Targeting the Anesthesia, Intensive Care, Emergency and Semi-Intensive Care markets, EDAN brings two new models iM60 and iM70 into the i-series family. These two new models are meant to add splendors to the i-Series portfolio and answer for the latest clinical requirements

Standard Accessories supplied with Basic Monitor.	Quantity
ECG Cable, 3 lead, snap, Defib, IEC, 3.5m, reusable	1
ECG Electrodes, adult, disposable, 30 pieces	30pcs
SPO ₂ Finger Sensor, adult, 2.5m, reusable (can be exchanged for child or Paed Probe)	1
NIBP Cuff, Adult	1
NIBP Cuff, Large Adult	1
NIBP Cuff, Child	1
NIBP Tube, 3m	1
Temperature Probe Skin, adult, 3m, reusable	1
Power Cord (3 Mtr)	1
Rechargeable Lithium-Ion Battery (14.8V, 4200mAh)	1
User Manual	1
Ground Cable	1

Multiple Display Modes



Standard Display



Large Font



Trend view



OxyCRG



Bed to Bed view

- 10.1" full-touch screen with optimized workflow
- Optional built-in Wi-Fi module
- Wired LAN with DHCP enabled
- HL7 support via XML files or MFM-CMS central monitoring system





iM70
Patient Monitor

EDAN's new generation patient monitor iM70 is suitable for both semi-intensive and intensive care. Thanks to compact design with touch screen, network capabilities, flexible parameters configuration as well as the reliable technology, iM70 will be the ideal choice for ER, OR, NICU, General Ward and also ideal for Patient Transport.

12.1"

RT2 Item Code
42181701-00001
Offer 3




Standard Parameters: 3/5-lead ECG, RESP, EDAN SpO2, EDAN NIBP, PR, 2-TEMP

Reliable Algorithms

- iSEAP™ ECG algorithm optimized for arrhythmia detection, pacemaker detection, and HR measurement
- iMAT™ SpO2 algorithm with outstanding motion resistance and low perfusion resistance performance
- iCUPS™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients

Multiple Display Modes



Standard Display Large Font Trend View OxyCRG Bed to Bed

Standard Accessories supplied with Basic Monitor.	Quantity
ECG Cable, 3 lead, snap, Defib, IEC, 3.5m, reusable	1
ECG Electrodes, adult, disposable, 30 pieces	30pcs
SPO2 Finger Sensor, adult, 2.5m, reusable (can be exchanged for child or Paed Probe)	1
NIBP Cuff, Adult	1
NIBP Cuff, Large Adult	1
NIBP Cuff, Child	1
NIBP Tube, 3m	1
Temperature Probe Skin, adult, 3m, reusable	1
Power Cord (3 Mtr)	1
Rechargeable Lithium-Ion Battery (14.8V, 4200mAh)	1
User Manual	1
Ground Cable	1


Full wall mounting and rolling stand options available

iM80
High-acuity Patient Monitor

The high performance of EDAN patient monitor iM80 is designed to match the needs for Critical Care, Post Anesthesia Care and Operation Theatres. The advanced monitoring system and built in clinical interfaces can provide the all-in-one solution for hospital segments

RT2 Item Code
42181701-00001
Offer 4

15"



Standard Parameters: 3/5-lead ECG, RESP, EDAN SpO2, EDAN NIBP, PR, 2-TEMP

- 15" high resolution colour display with maximum 13 waveforms
- Full touch-screen
- Pacemaker detection, electrosurgical interference proof
- SD card slot enables memory extension for large data storage
- USB and serial ports enables abundant future upgrades
- OxyCRG function for neonates
- Large font display
- Arrhythmia and S-T segment analysis
- Built-in rechargeable Li-ion Battery (up to 6 hours)
- Bi-directional communication with the MFM-CMS central station

Standard Accessories supplied with Basic Monitor.	Quantity
ECG Cable, 3 lead, snap, Defib, IEC, 3.5m, reusable	1
ECG Electrodes, adult, disposable, 30 pieces	30pcs
SPO2 Finger Sensor, adult, 2.5m, reusable (can be exchanged for child or Paed Probe)	1
NIBP Cuff, Adult	1
NIBP Cuff, Large Adult	1
NIBP Cuff, Child	1
NIBP Tube, 3m	1
Temperature Probe Skin, adult, 3m, reusable	1
Power Cord (3 Mtr)	1
Rechargeable Lithium-Ion Battery (14.8V, 4200mAh)	1
User Manual	1
Ground Cable	1

For comprehensive brochures, including full specifications, please see our website www.eastcoastmedical.co.za

Western Cape: 021 905 7721, Eastern Cape: 043 736 3111, KZN: 031 201 5585, Gauteng: 012 661 3513

**Monitor
Accessories
available on RT-2
Contract**



Pic	RT2 Item Code	Description
1	42181716-00001	CABLE, ECG, 3-LEAD, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
2	42181716-00002	PROBE, TEMPERATURE, WITH CABLE, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
3	42181716-00003	PROBE, SPO2, FINGER, WITH CABLE, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
4	42181716-00004	CUFF, NIBP, ADULT, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
5	42181716-00005	HOSE, NIBP, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
6	42181716-00006	LEAD, ECG, ADULT, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
7	42181716-00007	LEAD, ECG PAEDIATRIC, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
8	42181716-00008	CUFF, NIBP, REUSABLE, PAEDIATRIC, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
9	42181716-00009	CUFF, NIBP, REUSABLE, ADULT, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
10	42181716-00010	CUFF, NIBP, REUSABLE, OBESE, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
11	42181716-00011	CUFF, NIBP, REUSABLE, THIGH, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
12	42181716-00012	PROBE, SPO2, FINGER, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
13	42181716-00013	PROBE, SPO2, EAR, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
14	42181716-00014	PROBE, SPO2, NEONATAL, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
15	42181716-00015	CABLE, SPO2, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
16	42181716-00016	PROBE, TEMPERATURE, SURFACE, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
17	42181716-00017	PROBE, TEMPERATURE, RECTAL, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
18	42181716-00018	BRACKET, WALL, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER
19	42181716-00019	ROLL STAND, REQUIRED FOR MONITOR, BASIC, MULTI-PARAMETER

 <p>1 CABLE, ECG, 3 LEAD</p>	 <p>2 PROBE, TEMPERATURE</p>	 <p>3 PROBE, SPO2 FINGER</p>	 <p>4 CUFF NIBP ADULT ~ 27cm-35cm</p>
 <p>5 HOSE NIBP , 3m</p>	 <p>6 LEAD, ECG, ADULT</p>	 <p>7 LEAD, ECG, PAEDIATRIC</p>	 <p>8 CUFF, PAEDIATRIC ~ 20cm-28cm</p>
 <p>9 CUFF NIBP ADULT ~ 27cm-35cm</p>	 <p>10 CUFF NIBP OBESE ~ 32-43cm</p>	 <p>11 CUFF NIBP THIGH ~ 46-66cm</p>	 <p>12 PROBE, SPO2 FINGER</p>
 <p>13 SPO2 SENSOR, EAR, ADULT</p>	 <p>14 SPO2 SENSOR, NEONATAL</p>	 <p>15 CABLE, SPO2 (Interface cable for disposable probes)</p>	 <p>16 PROBE, TEMPERATURE, SURFACE</p>
 <p>17 PROBE, TEMPERATURE, RECTAL</p>	 <p>18 Bracket, Wall (Includes basket)</p>	 <p>19 ROLL STAND (Includes Monitor mounting plate)</p>	<p>For additional accessories, please contact your nearest East Coast Medical branch</p>



RT2 Item Code

42181701-00002

elite V6

Modular Patient Monitor



MONITOR, ICU, MODULAR

Product Introduction

Engineered specially for high-acuity divisions, the elite V6 is dedicated to bringing high-quality health care to intensive care and anesthesia monitoring, integrating world-leading technologies into one unit.

XM Module (supplies as standard, as per RT2 specification)
Standard: 3/5-lead ECG, NIBP, SpO₂ with Signal Intensity (S1), 2-TEMP, 12-lead ECG, 2-IBP.

Optional Modules:

- V-SpO₂ Module (Neillcor OxiMax™ SpO₂ with SatSeconds™)
- V-NIBP Module (Omron® NIBP) •
- V-IBP Module (Maximum 8-IBP with waveform overlapping function)
- V-C.O. Module (Thermodilution Cardiac Output)
- V-ICG Module (Impedance Cardiography)
- V-CO₂ Module (Respironics Mainstream/Sidestream, G2 Sidestream)
- V-RM Module (Respironics Respiration Mechanics)
- V-AG Module (Masimo Mainstream/Sidestream Anaesthetic Gas)
- V-BIS Module (Bispectral index)
- V-AG Module (Dräger Sidestream AG, Single Channel, without AAID, with O₂)



The following items are supplied as standard - included in the price for the **EDAN elite V6 Modular Monitor** (price as shown on the RT2 Pricing Schedule)

PART NO.	Item	Standard	Quantity
83.60.352009	V6 Main Unit	15" Color TFT Touchscreen	1
Standard accessories:			
PART NO.	DESCRIPTION		Quantity
01.13.36014	Power Cable (RSA - 3Mtr Standard)		1
01.13.114214	Ground Cable		1
	Check ERP		1
XM Module			
Item	Configurations	Quantity	
Including	XM Module 3/5-lead ECG, RESP,NIBP,2-TEMP,PR, EDAN SpO ₂ (with standard accessories)	1	
Including	12-lead ECG option Built-in 3/5/12-lead ECG module with standard accessories	1	
Including	2-IBP option Built-in 2-IBP module with standard accessories	1	
Standard accessories supplied with monitor:			
PART NO.	DESCRIPTION		Quantity
ECG			
01.57.471228013	3/5-lead ECG trunk Cable with Leadwires kit, Snap (IEC): 1 x ECG trunk cable, 5-lead, Defib, IEC, 2.6m, reusable		1
01.13.036627013	1 x ECG limb wires, 5-lead, snap, IEC, 0.9m, reusable		
01.57.471168014	12-lead ECG trunk Cable with Leadwires kit, Snap (IEC): 1 x ECG trunk cable, 10-lead, Defib, IEC, 2.6m, reusable		1
01.57.040203014	1 x ECG limb wires, 10-lead, snap, IEC, 0.9m, reusable		
01.57.471276-10	ECG Electrodes, adult, disposable, 10 pcs		10 pcs
SpO₂			
02.01.210122014	EDAN SpO₂ Accessories kit: 1 x SpO ₂ Silicone Soft-tip Sensor, adult, 1m, reusable		1
01.57.471068013	1 x SpO ₂ 7-pin Extension Cable		
NIBP			
01.57.471330011	EDAN NIBP Accessories kit: 1 x NIBP Cuff, Adult, 27cm-35cm, reusable, with EDAN logo		1
01.59.473007011	1 x NIBP Tube		
IBP			
01.57.471070014	Pressure transducer interface cable, BD		1
01.57.40121	IBP Pressure transducer kit, BD, disposable (BD DT-4812)		1
Temp			
01.15.040225012	2-pin Adult Skin Temperature Probe (10K/25°C)		1

ELECTROSURGICAL UNITS

RT2 Item Code

42295104-00025

ELECTROSURGERY UNIT, GENERAL SURGERY (350 W), SOLID STATE

OBS-350A



FDA 510(K)



Features:

- **Automatic set-up** at switch on, with fault detection
- **Memory function** for 10 users. Unit will start-up at settings when previously switched off
- **Patient Plate contact quality** inspection system - the contact quality between the patient and patient plate is measured in real time, ensuring continual patient safety.
- **Automatic activation** of Argon System, if matched up with an Argon Beam system.
- **User friendly**, with large displays for easy viewing
- **Excellent compatibility** with Cystoscopes, Arthroscopes, laparoscopes, etc - providing safe electrosurgery. Unit also ideal for TURP and TUVF procedures.
- **Independent output** terminals for each channel

Working Modes :

- Pure cut: 1w ~ 350w (Load 500Ω)
- Blend1: 1w ~ 250w (Load 500Ω)
- Blend2: 1w ~ 200w (Load 500Ω)
- Coagulation1: 1w ~ 150w (Load 500Ω)
- Coagulation2: 1w ~ 120w (Load 500Ω)
- Bipolar: 1w ~ 70w (Load 100Ω)

Standard accessories:

Disposable electrosurgical pencil, disposable pad, plate cable, reusable bipolar forceps, footswitch, mains cable.

Dimensions: 47cm*42cm*16cm

Gross weight: 13Kg

RT2 Item Code

42295104-00039

EVACUATION, SMOKE, FOR ELECTROSURGERY UNIT, GENERAL SURGERY



Pic	RT2 Item Code	Accessories
1	42295104-00026	FOOTSWITCH, DOUBLE FOR ELECTROSURGERY UNIT, GENERAL SURGERY
2	42295104-00027	ELECTRODE, REUSABLE STAINLESS STEEL OR CONDUCTIVE RUBBER BIPOLAR PATIENT ELECTRODE FOR ELECTROSURGERY UNIT, GENERAL SURGERY
3	42295104-00028	CABLE, PATIENT ELECTRODE FOR ELECTROSURGERY UNIT, GENERAL SURGERY
4	42295104-00029	HANDLES, REUSABLE, STERILIZABLE, ELECTRODE WITH CABLES, 2 SETS FOR ELECTROSURGERY UNIT,
5	42295104-00030	ELECTRODES, BLADE/KNIFE, REUSABLE, STRAIGHT FOR ELECTROSURGERY UNIT, GENERAL SURGERY
6	42295104-00031	ELECTRODES, NEEDLE, THICK, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL SURGERY
7	42295104-00032	ELECTRODES, NEEDLE, THIN FOR ELECTROSURGERY UNIT, GENERAL SURGERY
8	42295104-00033	ELECTRODES, BALL, STRAIGHT, 4 MM DIAMETER, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL SURGERY
9	42295104-00034	ELECTRODES, BALL, STRAIGHT, 2,5 MM DIAMETER FOR ELECTROSURGERY UNIT, GENERAL SURGERY
10	42295104-00035	ELECTRODES, BIPOLAR, REUSABLE, STERILIZABLE, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL
11	42295104-00036	PAD, ELECTRODES, DISPOSABLE SINGLE OR SPLIT NEUTRAL FOR ELECTROSURGERY UNIT, GENERAL
12	42295104-00037	SEALING, VESSEL, FOR ELECTROSURGERY UNIT, GENERAL SURGERY





Electrosurgical Accessories

available on RT2 Contract

CE ISO 9001 13485 FDA 510(K)

1



RT2 Item Code:
42295104-00026

FOOTSWITCH, DOUBLE FOR ELECTROSURGERY UNIT, GENERAL SURGERY

2



RT2 Item Code:
42295104-00027

CONDUCTIVE RUBBER PATIENT PLATE - REUSABLE

3



RT2 Item Code:
42295104-00028

CABLE, PATIENT ELECTRODE

5



RT2 Item Code:
42295104-00030

ELECTRODES, BLADE/KNIFE, REUSABLE, STRAIGHT FOR ELECTROSURGERY UNIT, GENERAL SURGERY

4 RT2 Item Code : 42295104-00029

HANDLES, REUSABLE, STERILIZABLE, ELECTRODE WITH CABLES, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL SURGERY




6



RT2 Item Code:
42295104-00031

ELECTRODES, NEEDLE, THICK, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL SURGERY

7



RT2 Item Code:
42295104-00032

ELECTRODES, NEEDLE, THIN FOR ELECTROSURGERY UNIT, GENERAL SURGERY

8



RT2 Item Code:
42295104-00033

ELECTRODES, BALL, STRAIGHT, 4 MM DIAMETER, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL SURGERY

9



RT2 Item Code:
42295104-00034

ELECTRODES, BALL, STRAIGHT, 2.5 MM DIAMETER FOR ELECTROSURGERY UNIT, GENERAL SURGERY

10



#0550-01 Length: 100mm Tip size: 0.5mm



#0550-05 Length: 110mm Tip size: 0.5mm



#0551-01 Length: 200mm Tip size: 0.5, 1.0, 1.5, 2.0mm



#0550-06 Length: 115mm Tip size: 0.5mm



#0550-08 Length: 195mm Tip size: 0.5, 1.0, 1.5, 2.0mm




#0551-05 Length: 200mm Tip size: 0.5, 1.0, 1.5, 2.0mm


RT2 Item Code: 42295104-00035

ELECTRODES, BIPOLAR, REUSABLE, STERILIZABLE, 2 SETS FOR ELECTROSURGERY UNIT, GENERAL SURGERY (User to select any 2 of the above Bipolar Electrodes)


11




#1030a Foam backing Adult, Mono-polar, Upright




#1030b Foam backing Child, Mono-polar, Upright




#1030n Foam backing Newborn, Mono-polar, Upright



#1031a Foam backing Adult, Bi-polar, Upright



#1031b Foam backing Child, Bi-polar, Upright



#1031n Foam backing Newborn, Mono-polar, Upright

RT2 Item Code: 42295104-00036

PAD, ELECTRODES, DISPOSABLE SINGLE OR SPLIT NEUTRAL FOR ELECTROSURGERY UNIT, GENERAL SURGERY. (User to select any one of the patient plates shown)

12



RT2 Item Code: 42295104-00036

SEALING, VESSEL, FOR ELECTROSURGERY UNIT, GENERAL SURGERY

Unit supplied complete with autoclaveable bipolar cable

RT2 Item Code

42191803-00011



YP-920

INCUBATOR ICU

- Air or Servo controlled Incubator;
- Humidity concentration control system;
- Set temperature, air temperature, baby temperature and humidity concentration are displayed separately;
- >37 °C temperature set function;
- Failure alarm indication;
- Embedded integrated sensor box;
- Integrated water reservoir;
- Double wall hood for temperature stabilisation;
- Inclination of bassinet is adjustable;
- RS-232 connector;
- Weighing system and Monitor shelf are optional.

RT2 Item Code

42191812-00001



YP-910

INCUBATOR (STANDARD)

- Air mode and baby mode controlled by micro-computer;
- Humidity concentration display system;
- Set temperature, air temperature, baby temperature, humidity concentration are displayed separately;
- >37°C temperature set function;
- Failure alarm indication;
- Embedded integrated sensor box;
- Double wall hood, automatic air circulation device;
- Inclination of bassinet is adjustable;
- RS-232 connector.
- Weighing system and Monitor shelf are optional

Item Code	Item Description
42191803-00012	FILTER, BACTERIAL, FOR INCUBATOR ICU
42191803-00013	SLEEVES, IRIS PORT, FOR INCUBATOR ICU, BOX OF 10
42191803-00014	SLEEVES, OVAL PORT, FOR INCUBATOR ICU, BOX OF 10
42191803-00015	PROBE, SKIN, TEMPERATURE, REUSABLE, FOR, INCUBATOR, ICU
42191803-00016	COVERS, TEMPERATURE, PROBE, FOR, INCUBATOR, ICU
42191803-00017	BATTERY, INTERNAL ALARMS, FOR INCUBATOR ICU
42191803-00018	SCALE, FOR INCUBATOR ICU
42191803-00019	LAMP, OBSERVATION, FOR INCUBATOR ICU

BRANCHES

Gauteng

Gateway Industrial Park
3 Display Gardens
Sarel Baard Crescent
Rooihuiskraal
Pretoria

Tel: 086 100 0428 / 012 661 3513
Fax: 012 661 7806
Cell: 082 573 9131
Email: sales@eastcoastmedical.co.za

KwaZulu Natal

East Coast Medical
305 Moore Road (Che Guevara)
Berea
Durban

Tel: 031 201 5585
Fax: 031 202 2442
031 201 1485
Cell: 082 485 0565
Email: james@eastcoastmedical.co.za

Cape Town

6 Pontac Road
Saxenburg Park 1
Blackheath
Cape Town

Tel: 021 905 7721
Fax: 021 905 7726
Cell: 083 265 8660
Email: renier@eastcoastmedical.co.za

East London

31 Sunningdale Avenue
Sunnyridge
East London

Tel: 043 736 3111
Cell: 083 272 9551
Fax: 043 736 3959
Email: lori@eastcoastmedical.co.za