

RT2-1-2016

Product Catalogue



**Bid Description: Supply, Delivery, Installation,
Commissioning and Maintenance of Medical
Equipment to the State**

INDEX

Category	Page
Blood Gas Analyzer i15	1
Tourniquet System	2
Colposcope	3
Echocardiogram AX4	4
Echocardiogram AX8	5



RT2-1-2016 BLOOD GAS ANALYZER EDAN i15 POCT

Line Item#	13.00
Item Code	41115806-00001
Item Description	BLOOD GAS ANALYZER Includes:100x BG10 Test Cartridges,1x Calibrant Pack, 10x Paper Rolls,1x UPS and Mobile Stand
Brand & Model	Edan i15 POCT
Lead Time	56 Days
Warranty	24 Months

Additional Accessories:

- Calibrant/Pack 1x CP100
- Test Cartridge- 100x BG10
- BGE 1 Quality Control
- BGE 3 Quality Control
- Electronic Simulator
- Printer Paper 1x Roll
- ECM UPS



Product Introduction: i15 Blood Gas and Chemistry Analyzer (POCT)

As a brand new POCT device specially designed for critical care, i15 is portable, easy to use and performs a flexible panel tests from a single disposable cartridge. i15 brings new era of blood gas and chemistry analysis, and it enables you to more quickly and effectively monitor and manage your patients.

Convenient operation, accurate and reliable results

- Standby without power calibrant wastage and maintenance free
- Fast, accurate, and convenient built-in auto calibration
- Auto sample aspiration assures convenience and reliability

Main features:

Small, fast, easy to use, cost efficient

BG10 Test Parameters

- Measured parameters: PH,PCO2,PO2, Na+,K+,Ca++,Cl-,Hct,Glu,Lac
- Calculated parameters: cH+, HCO3-act, HCO3-std, BE(ecf), BE(B), BB(B), ctCO2, sO2(est),
- Calculated parameters: Ca++(7.4), AnGap, tHb(est),pO2(A-a), pO2(a/A),RI, pO2/FIO2, cH+(T), pH(T),
- Calculated parameters: pCO2(T), pO2(T), pO2(A-a)(T), pO2(a/A)(T), RI(T), pO2(T)/FIO2 , mOsm

RT2-1-2016

TOURNIQUET SYSTEM DS MAREF DTS-3000

Line Item#

41.00

Item Code

42295122-00001

Item Description

Tourniquet System Includes: Single Cuffs (1x Paeds, 1x Adult), Double Cuffs(1x Paeds, 1x Adult,1x Large) and Mobile Stand with Basket.

Brand & Model

DS MAREF DTS-3000

Lead Time

56 Days

Warranty

24 Months

Additional Accessories:

Tourniquet Cuff Single 40 x 7cm

Tourniquet Cuff Single 61 x 9cm

Tourniquet Cuff Single 86 x 10cm

Tourniquet Cuff Double 57 x 10cm

Tourniquet Cuff Double 80 x 15cm

Tourniquet Cuff Double 86 x 15cm

Bandage, Reusable, Esmarch

Tourniquet Cuff Main Hose 200cm

Cuff Hose 15cm



DTS-3000w systems...

- 2 channels of individual settings & operations
- 2 channels (4 Ports for 2 double or single cuffs use)
- Easy settings by Touch Key operation
- LCD screen : Easy checking and cleaning
- Save function
- Operation recording function by SD memory card
- LOP sensor for automatic optimum pressure setting
- Self air pressure sensor calibration
- Safety alarm system (cuff, pressure, battery, time)
- Compressor dimension: 225(W) X 185(D) X 240 (H)mm
- Built-in battery (Max life time: 12 ~72 hours)
- Max pressure: 600mmHg
- Voltage: Free voltage (AC 100 - 240V~ 50/60Hz)
- Power consumption: 50W
- Noise dB: Below 60dB
- Max pressure: 600mmHg
- Time range: 5~240 Min.
- Compressor weight: 5.6kg
- Self test of pressure and system check

RT2-1-2016

COLPOSCOPE EDAN C6 HD

Line Item# 49.00
Item Code 42296309-0000


Item Description **COLPOSCOPE** Includes: C6 HD Screen, Swing Arm Trolley, Wireless Mouse, Dust Covers, Printer with 1000x sheets paper, 1x Colour Ink Cartridge and R-Way Cervical Screening Software Package.

Brand & Model EDAN C6 HD
Lead Time 56 Days
Warranty 24 Months

Additional Accessories:
 Colposcope System Dust Cover
 Wireless Mouse
 Printer
 Keyboard
 A4 or A6 Paper (1000 sheets)
 PC Software
 Camera Dust Cover



R-way™ is a patented evaluation program which assists in auto diagnosis according to preset standardized Colposcopic screening procedures to quickly identify CIN2+ cervical diseases.



RT2-1-2016

ECHOCARDIOGRAM MACHINE - PORTABLE EDAN ACCLARIX AX4

Line Item#

63.00

Item Code

42201709-00005

Item Description

Echocardiogram Machine - Portable Includes: 1x P5-1Q adult echo transducer, Battery Backup and 1x Multi-Transducer Connector

Brand & Model

Edan Acclarix AX4

Lead Time

56 Days

Warranty

24 Months

Additional Accessories:

Adult Echo Transducer

Transthoracic

Echocardiography Transducer

Innovative Linear

3D/4D Curved Transducer

External DVD/CD Writer

Paeds Echo Transducer

Vaginal Probe

Obstetric Probe

B/W Video Printer

MT-807 Trolley

Travel Case



System Description:

The remarkable Acclarix AX4 Compact Ultrasound System delivers a power house combination of features to meet the demands of point-of-care and general imaging applications. The Acclarix AX4 has been designed from the ground up with relentless focus on delivering unexpected levels of innovation and performance at a price point that is equally surprising.

- * Sleek compact design for portable exam environments.
- * Customizable control panel to personalize your user interface.
- * Sealed control panel to aid with infection control.
- * Tilt-and-swivel monitor for optimal viewing during procedures.

RT2-1-2016

ECHOCARDIOGRAM MACHINE - PORTABLE EDAN ACCLARIX AX8

Line Item#	63.00
Item Code	422017109-00005
Item Description	<i>Echocardiogram Machine - Portable</i> Includes: 1x P5-1Q adult echo transducer, Battery Backup and 1x Multi-Transducer Connector
Brand & Model	Edan Acclarix AX8
Lead Time	56 Days
Warranty	24 Months

Additional Accessories:

- Adult Echo Transducer
- Transthoracic Echocardiography Transducer
- Innovative Linear
- 3D/4D Curved Transducer
- External DVD/CD Writer
- Paeds Echo Transducer
- Vaginal Probe
- Obstetric Probe
- B/W Video Printer
- MT-807 Trolley
- Travel Case



Acclarix AX8 -- New Generation Ultrasound Platform- Silicon Valley innovation and design expertise 128 channels on transmit and receive Exceptional computing power resulting in outstanding image processing The Acclarix AX8 System boasts a full range of transducers that span most point-of-care and general imaging applications. All transducers feature multiple, selectable fundamental and harmonic frequencies to provide exquisite sensitivity, excellent detail and contrast resolution, and image uniformity throughout the field of view so you can appreciate even the smallest details.

- * Sleek compact design for portable exam environments.
- * Customizable control panel to personalize your user interface.
- * Sealed control panel to aid with infection control.
- * Tilt-and-swivel monitor for optimal viewing during procedures.

BRANCHES

www.eastcoastmedical.co.za

GAUTENG

Gateway Industrial Park
3 Display Gardens
Sarel Baard Crescent
Rooihuiskraal

Tel: 086 100 0428 / 012 661 3513
Fax: 012 661 7806
Cell: 082 573 9131
Email: sales@eastcoastmedical.co.za

KWAZULU NATAL

East Coast Medical
305 Moore Road (Che Guevara)
Berea
Durban

Tel: 031 201 5585
Fax: 031 202 2442
Fax: 031 201 1485
Email: james@eastcoastmedical.co.za

CAPE TOWN

6 Pontac Road
Saxenburg Park 1
Blackheath
Cape Town

Tel: 021 905 7721
Fax: 021 905 7726
Cell: 083 265 8660
Email: renier@eastcoastmedical.co.za

EAST LONDON

31 Sunningdale Avenue
Sunnyridge
East London

Tel: 043 736 3111
Fax: 043 736 3959
Cell: 083 272 9551
Email: lori@eastcoastmedical.co.za