

# East Coast Medical

**Northern Regions**

## RT2-2021 CATALOGUE

NATIONAL TREASURY TRANSVERSAL CONTRACT ENDING NOVEMBER 2024

Northern Regions - Pretoria  
Tel: 086 100 0428, 012 661 3513/7565/893  
Email: [sales@eastcoastmedical.co.za](mailto:sales@eastcoastmedical.co.za)

Kwazulu Natal - Durbn  
031 201 585  
Email: [james@eastcoastmedical.co.za](mailto:james@eastcoastmedical.co.za)



[www.eastcoastmedical.co.za](http://www.eastcoastmedical.co.za)

ITEM DESCRIPTION**Intubation Video Disposable Flexible Bronchoscope**

Item Image

Model

**SEESHEEN****Main Offer: BR-1231 Video Broncoscope System**

Offer includes:

- ▶ EC5311 Cat with IV pole
- ▶ 8" control panel. PAL compatible
- ▶ 1x flexible bronchoscope

**Category 2: RT2-02-09**ITEM DESCRIPTION**Resuscitator, Gas-powered, Infant**

Item Image

Model


**Offer 1: Neo-1**

Offer includes:

- ▶ Starter pack of consumables

**Category 2: RT2-02-11**ITEM DESCRIPTION**Air-O2 Blender, Low Flow**

Item Image

Model


**Main Offer: KY-20A**

Offer includes:

- ▶ Colour coded air & oxygen hoses with Heyer-type connectors
- ▶ Mounting bracket
- ▶ Flow meter with control knobs



ITEM DESCRIPTION**ICU Crib**

Item Image

Model

**Main Offer: HKN-93C**

- Offer includes:
- ▶ Skin temperature probe (reusable)
  - ▶ Disposable hydrophobic, anti-bacterial filters. Box of 10
  - ▶ Air / oxygen blender
  - ▶ LED phototherapy
  - ▶ IV pole
  - ▶ Shelf for monitor
  - ▶ Instrument tray - Internal
  - ▶ Internal battery for alarms

**Category 3: RT2-03-08**ITEM DESCRIPTION**Incubator, ICU**

Item Image

Model

**Offer 1: YP-920**

- Offer includes:
- ▶ Skin temperature probe (reusable)
  - ▶ Covers for temperature probe
  - ▶ Anti-bacterial filters x4
  - ▶ Internal battery for alarms

**Category 3: RT2-03-09**ITEM DESCRIPTION**Incubator, Standard**

Item Image

Model

**Offer 2: YP-910**

- Offer includes:
- ▶ Skin temperature probe (reusable)
  - ▶ Covers for temperature probe x25
  - ▶ Anti-bacterial filters x4
  - ▶ Internal battery for alarms
  - ▶ Oval port sleeves x10
  - ▶ Iris port sleeves x10

ITEM DESCRIPTION**Incubator, Transport**

Item Image

Model

**Main Offer: TI-2000**

Offer includes:

- ▶ Skin temperature probe (reusable)
- ▶ Covers for temperature probe x25
- ▶ Anti-bacterial filters x4

**Category 3: RT2-03-11**ITEM DESCRIPTION**Infant Resuscitator**

Item Image

Model

**Main Offer: HKN-93B**

Offer includes:

- ▶ Skin temperature probe (reusable)
- ▶ Mattress
- ▶ IV pole
- ▶ Accessory rail
- ▶ Monitor stand x2

**Category 3: RT2-03-12**ITEM DESCRIPTION**Phototherapy Unit, LED**

Item Image

Model

**Main Offer: XHZ-90P**

Offer includes:

- ▶ Phototherapy unit
- ▶ All inclusive LED examination light attached to mobile stand



## ITEM DESCRIPTION

## Ultrasound Machine - Neonatal

Item Image

Model



## Offer 2: Acclarix LX9

Offer includes:

- ▶ Mobile cart based
- ▶ Linear array transducer: L12-5Q
- ▶ Phased array transducer: P5-1Q
- ▶ Phased array transducer: P7-3Q
- ▶ Battery (5000mAh)
- ▶ Ultrasound gel
- ▶ Sony UP-X898MD BW printer
- ▶ Sony UP-X898MD BW printer paper 10 rolls/pack
- ▶ 4Hr UPS

## Category 4: RT2-04-03

## ITEM DESCRIPTION

## Cardiotocograph, Twins on Trolley

Item Image

Model



## Offer 3: FTS-6

Offer includes:

- ▶ FTS-6 FHR probe (Ultrasound probe) x8
- ▶ FTS-6 TOCO probe x7
- ▶ FTS-6 Belt x15
- ▶ Rechargeable lithium-ion battery
- ▶ Ultrasound gel
- ▶ ECM colour printer + ink cartridges + pack of A4 paper x2
- ▶ Wireless mouse
- ▶ 4Hr UPS
- ▶ 3m kettle plug
- ▶ Mobile cart based

## Category 4: RT2-04-04

## ITEM DESCRIPTION

## Cardiotocograph, Twins + MHR on Trolley

Item Image

Model



## Offer 1: F15 twins, TOCO-MHR (SPO2, ECG, NIBP)

Offer includes:

- ▶ FHR probe (Ultrasound probe) x2
- ▶ F15 TOCO probe / SpO2 probe / ECG cable / NIBP hose / Event marker
- ▶ F15 Belt x3
- ▶ F15 Ultrasound gel
- ▶ F15 Recording paper x5
- ▶ F15 Rechargeable lithium-ion battery
- ▶ Trolley
- ▶ 3m kettle plug



ITEM DESCRIPTION**Fetal Doppler - Deep Penetration For Later Pregnancy**

Item Image

Model

**Main Offer: Edan SD3 Pro, No Trolley**

Offer includes:

- ▶ 2Mhz waterproof interchangeable probe
- ▶ Rechargeable batteries
- ▶ Charging stand
- ▶ Carry case

**Category 4: RT2-04-06**ITEM DESCRIPTION**Video Colposcope and Accessories**

Item Image

Model

**Main Offer: C6 HD on Trolley with Swing Arm**

Offer includes:

- ▶ Colour printer
- ▶ Colour ink cartridges: colours x3, black x1
- ▶ Pack of A4 paper x2
- ▶ Wireless mouse
- ▶ Dust cover
- ▶ 3m kettle plug
- ▶ C6 swing arm stand

**Category 5: RT2-05-07**ITEM DESCRIPTION**ECG Machine, 12-Lead**

Item Image

Model

**Offer 1: SE-1200 Express with accessories and trolley**

Offer includes:

- ▶ MT-801 trolley
- ▶ 4mm banana ECG cable with crocodile clips
- ▶ Rechargeable lithium-ion battery (1600mAH)
- ▶ Adult long term electrodes 50pcs/pack

ITEM DESCRIPTION**Patient Monitor, Basic**

Item Image

Model

**Offer 1: iM80 (3/5-lead ECG, RESP, NIBP, SpO2, 2-TEMP)**

Offer includes:

- ▶ Paeds 3/5 lead ECG cable
- ▶ Adult 3/5 lead ECG cable
- ▶ ECG interface cable
- ▶ Reusable latex-free, neonatal, paedes, adult, obese cuff
- ▶ NIBP hose
- ▶ SpO2 finger sensor for paedes & adult
- ▶ SpO2 wrap probe for neonatals
- ▶ SpO2 interface cable
- ▶ 2x reusable temperature probes

**Category 5: RT2-05-15**ITEM DESCRIPTION**Patient Monitor, ICU**

Item Image

Model

**Offer 1: V8 (3/5 lead ECG, RESP, NIBP, SpO2, 2-TEMP with XM module with IBP)**

Offer includes:

- ▶ Paeds 3/5 lead ECG cable
- ▶ Adult 3/5 lead ECG cable
- ▶ ECG interface cable
- ▶ Reusable latex-free, neonatal, paedes, adult, obese cuff
- ▶ NIBP hose
- ▶ SpO2 finger sensor for neonatals, paedes & adult
- ▶ SpO2 interface cable
- ▶ Reusable temperature probes x2
- ▶ Offer upgradable with CO2

**Category 5: RT2-05-16**ITEM DESCRIPTION**ICU Patient Monitor with Transport Module**

Item Image

Model

**Offer 1: V8 ICU monitor with iM20 module (3/5 Lead ECG, RESP, NIBP, SpO2, 2-TEMP, Incl.: 2-IBP & 12-Lead ECG)**

Offer includes:

- ▶ Paeds & Adult 3/5 lead ECG cable
- ▶ ECG interface cable
- ▶ Reusable latex-free, neonatal, paedes, adult, obese cuff
- ▶ NIBP hose
- ▶ SpO2 finger sensor for paedes & adult
- ▶ SpO2 wrap sensor for neonatals
- ▶ SpO2 interface cable
- ▶ Reusable temperature probes x2

ITEM DESCRIPTION**Vital Signs Monitor with Oximetry & Quick Temp**

Item Image

Model

**Offer 1: iM3-NST with trolley, includes all adult accessories**

Offer includes:

- ▶ Reusable latex-free paed's / neonatal / adult / obese cuff
- ▶ NIBP hose
- ▶ Finger sensor for paed's / adults
- ▶ SpO2 interface cable
- ▶ SpO2 wrap probe for neonatals
- ▶ Reusable temperature probes x2
- ▶ MT-206 trolley

**Category 5: RT2-05-19**ITEM DESCRIPTION**Vital Signs Monitor with Oximetry & Quick Temp**

Item Image

Model

**Offer 2: M3A-NST with trolley, includes all adult accessories**

Offer includes:

- ▶ Reusable latex-free paed's / neonatal / adult / obese cuff
- ▶ NIBP hose
- ▶ Finger sensor for paed's / adults
- ▶ SpO2 interface cable
- ▶ SpO2 wrap probe for neonatals
- ▶ Reusable temperature probes x2
- ▶ MT-206 trolley

**Category 7: RT2-07-05**ITEM DESCRIPTION**Tourniquet System**

Item Image

Model

**Main Offer: DTS-3000**

Offer includes

- ▶ Single cuff: paed's, adult + double cuff: paed's, adult & large
- ▶ All hoses
- ▶ Tourniquet cuff main hose: 200cm
- ▶ Cuff hose 15cm
- ▶ Reusable, Esmarch bandages
- ▶ Rolling stand with basket
- ▶ Integrated rechargeable battery backup
- ▶ LOP sensor for automatic, optimum pressure settings



## ITEM DESCRIPTION

## Laparoscope Instruments

Item Image

Model

## Offer: SonoScape Accessories



Straight Forward Telescope 0° enlarged view, diameter 10mm, length 29-31cm, autoclavable, fiber optic light transmission incorporated	
Straight Forward Telescope 0° enlarged view, diameter 5mm, length 29-31cm, autoclavable, fiber optic light transmission incorporated	
Straight Forward Telescope 0° enlarged view, diameter 2.5mm, length 29-31cm, autoclavable, fiber optic light transmission incorporated	
Forward-Oblique Telescope 30° enlarged view, diameter 10mm, length 29-31cm, autoclavable, fiber optic light transmission incorporated	
Forward-Oblique Telescope 30° enlarged view, diameter 5mm, length 29-31cm, autoclavable, fiber optic light transmission incorporated	
f Forward-Oblique Telescope 30° enlarged view, diameter 2.5mm, length 29-31cm, autoclavable, fiber optic light transmission incorporated	
g Plastic Container for Sterilizing and Storage, perforated, with lid for two-level storage, for use with 30cm and 36cm instruments, external dimensions (w x d x h): 600mm x 145mm x 260mm	
h Laparoscopic Needle Holder	
i VERESS Pneumoperitoneum needle, with spring-loaded blunt stylet, LUER-Lock, length 10cm	
j Trocar size, 11mm, with insufflation stopcock, Cannula with thread, length 100mm (max 110mm)	
k Trocar size, 3.5mm, with insufflation stopcock, Cannula with thread, length 100mm (max 110mm)	
l Trocar size, 6mm, with insufflation stopcock, Cannula with thread, length 100mm (max 110mm)	

NB: \*Prices indicated herein, are per unit &amp; include VAT. Subject to annual approved escalation.

Effective upto 30 November 2024